MINIPIX

Version 1.0

Datasheet

General description

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
Illustration of single particle sensitivity of Timepix device. The tracks of different particles of normal radiation background were recorded in 10minute. Nonoise (clean zero) is seen in dark regions.

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¹ MINIPIX is not certified dosimetric device. It serves as the first level indicator and monitor of radiation fields allowing identification of a radiation type. Radiation protection of people cannot be based on measurements of MINIPIX.

² Dynamic range of final picture is theoretically unlimited; the only limiting factor is exposure time.

Device Parameters

Operating Conditions

Symbol	Parameter	Min	Typ	Max	Units	Comment
T_A	Temperature Range	0	50	70	°C	
Φ	Humidity			60	%	Not condensing

Electrical Specification

$T_A = 25^\circ\text{C}$, USB voltage $V_{CC} = 4.8\text{V}$


Symbol	Parameter	Min	Typ	Max	Units	Comment
V_{CC}	Supply Voltage	4.0	5.0	6.0	V	Comply with USB 2.0
I_{CC}	Supply Current					
I_{CC1}	Chip disabled		147.0		mA	
I_{CC2}	Chip active			500	mA	Comply with USB 2.0
P_2	Power Dissipation			2.5	W	
V_{INL}	Aux. Conn. Input Voltage Low	-0.3		1.16	V	
V_{INH}	Aux. Conn. Input Voltage High	1.83		5.5	V	
Bias Voltage Source for Sensor Diode						
V_{BIAS}	Bias Voltage	3		200	V	
I_{BIAS}	Bias Current ($V_{BIAS} = 5\text{V}$)			10	mA	
	Bias Current ($V_{BIAS} = 30\text{V}$)			5	mA	
	Bias Current ($V_{BIAS} = 100\text{V}$)			1	mA	

Performance Characteristics

Symbol	Parameter	Min	Typ	Max	Units	Comment
f	Frame Rate			30	fps	with USB 2.0 Host
T_{READ}	Frame Readout Time ³		33		ms	

³ During Readout time (or Dead time), no charge is collected from the sensor.

Device description

USB connector





USB type mini B, Standard USB 2.0 Hi-Speed.

Auxiliary connector

Connector is available internally on board. It can be made accessible upon request.

1	GND	2	+5V
3	Shutter	4	Reserved
5	Clock	6	UART TX
7	Reserved	8	UART RX
9	Reserved	10	Reserved

Mechanical Dimensions

All dimensions are in mm

Warning

- **Do not touch sensor surface!**

Instructions for safe use

To avoid malfunction or damage to your **MINIPIX** please observe the following

- Do **not** expose to water, moisture.
- Do **not** disassemble. Wire-bonding connection may be irreversibly damaged.
- Do **not** insert any object into the sensor window.

Copyright

WIDEPIX s.r.o.

Na Balkane 2075/70, CZ 130 04 Praha 3

Czech Republic

Tel: +420-603-444112, 589854;

Email: info@widepix.cz

www.widepix.cz

For more information visit **WIDEPIX**[®] website at www.widepix.cz.

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